

Sixteen Treatment Protocols to Choose From

The NM III™ is a versatile, fully adjustable muscle stimulator that is capable of delivering two unique waveforms, has 16 preset programs and can be manually programmed. This device can be tailored to suit the numerous and individual needs of patients needing muscle strengthening and gains in range of motion.

THE NM III OFFERS THESE OUTSTANDING BENEFITS

Choice of Two Waveforms

The modified rectangular asymmetric waveform option on the NM III is very effective in treating conditions involving smaller muscles and tendons of the extremities. It allows the clinician to isolate and treat a specific muscle or small muscle groups. This isolation minimizes fatigue and eliminates unwanted muscle responses.

The symmetrical biphasic waveform comfortably delivers the stimulation necessary for the larger muscle groups of the upper legs, trunk and torso. This waveform stimulates the larger muscle groups, thus enhancing range of motion more quickly than volitional movements alone. The stimulation is equal under each electrode, improving comfort and enhancing patient compliance.

Choice of 16 preset programs

- Treats a wide range of conditions from ACL repair to muscle spasticity
- Easy to use
- Ideal for those who appreciate traditional protocol parameters

Flexible manual set-up options

- Lag-time selection helps restore normal muscle tracking and timing
- Unlike other devices, the NM III features hand and heel switches which can be utilized in one or two channels
- Ideal for those who need more flexibility in their protocol design

Portable

The NM III is portable, wearable and easy-to-use. It is an ideal device for enhancing muscle rehabilitation in the home, the clinic or in the field.

COMMON USES

- Post ACL repair
- Muscle spasms
- Total Knee Replacement
- Shoulder Subluxation



III™

 **Rehabiliticare**

